## **DNA Ladders**

Omega Biotek offers a wide range of products for the accurate size and mass determination of DNA. Most of DNA Ladder includes one or two higher-intensity bands for easy identification and orientation. The ladders are supplied premixed with loading buffer and are stable at room temperature for a few weeks. There is no need to heat, mix, or dilute ladders prior to loading them on gel.

## 50bp DNA Ladder

CAT.No.: M11-01Packing: 250 ulConcentration:

• Store: at  $4^{\circ}$ C for periods up to 6 months. For longer periods store at  $-20^{\circ}$ C.

• **Description:** 50 DNA Ladder is designed for sizing and approximate quantification of wide range linear double-stranded DNA fragments on agarose gel. The ladder is composed of eight purified individual DNA fragments (in base pairs):500,400,300,250,200,150,100,50. It contains one reference band (250 bp) for easy identification and orientation.

## • Direction:

- ♦ This product can load directly, since it contains 1×loading buffer. It is free of fluorescence dye such as EtBr.
- ♦ For a loading well less than 6mm in width, 5ul of this product is enough for clear bands. For a wider well, more product may be required.
- ♦ Gel and buffer may effect the electrophoresis. For better result, please use high quality agarose and fresh running buffer. 3.0% agarose gel is recommended.

## Electrophoresis Map:

